

# Direct Monitoring of Combinatorial Chemistry Reactions by Infrared Microspectroscopy

Pamela A. Martoglio, Ph.D. and David W. Schiering, Spectra-Tech, Inc., Matthew J. Smith, Ph.D. Nicolet Spectroscopy Research Center; and Daniel T. Smith, The Pennsylvania State University.

**What is combinatorial chemistry?** Combinatorial chemistry is quickly becoming a very popular technique for organic synthesis in the pharmaceutical and biotechnology fields. Basically, combinatorial methods allow one to obtain thousands of derivatives of chemical compounds in a very quick and efficient manner. The method begins when a molecule of interest is attached to resin beads. Once the molecule is attached, a number of reactions can be run on the beads. The beads can be split up into several subsets, and different reactions can be run on each of the subsets. Because the beads are in the solid phase, it is extremely easy to separate the beads from the liquid layer after the reactions have completed.

After the first set of reactions is finished, the subsets can be mixed together and re-divided into a second group of subsets. Now, each subset contains several different reaction products of the molecules attached to the beads. By running different reactions on each of these subsets, one obtains even more reaction products. If this cycle is repeated, it becomes obvious how literally thousands of different derivatives of a single molecule can be combined rather quickly. Hence, this process earned the name of Combinatorial Chemistry.

**How can the reaction products be detected and characterized?** There are a variety of methods currently being investigated for the analysis of the final reaction products. Some groups cleave the molecule from the bead after the reaction, and detect it by traditional methods. We propose that infrared microspectroscopy could be used to detect and characterize the reaction products of the molecules, directly on a single bead, without having to cleave the molecule from the bead. This direct sampling capability will greatly speed the entire process of combinatorial chemistry.




Figure 1: Video image capture of a single unreacted combinatorial bead

In this paper, we present the results from three different bead samples: an unreacted bead (Figure 1) and two reacted beads. We will show that high-quality infrared spectra can be obtained of single beads, and the reaction products can easily be detected and characterized.


## Experimental

The unreacted beads (from Nova Bio Chem) were polystyrene based, with a polyethylene glycol group added on to the polystyrene base. Attached to the polyethylene glycol group was an active amine group. This amine group was the reactive site where core molecules were added. The reactions are detailed below.



Scheme 1

To a mixture of 4.86 mg of beads in  $\text{CH}_2\text{Cl}_2$  (1 ml) was added 50  $\mu\text{l}$  of pyridine (M.W. = 79.10) and 50  $\mu\text{l}$  of acetic anhydride (M.W. = 102.19). The mixture was stirred at room temperature ( $25^\circ\text{C}$ ) for 24 hours. The mixture was filtered (through Whatman qualitative filter paper cat. 1001 042), washing five times with  $\text{CH}_2\text{Cl}_2$ . The beads and filter paper were put under vacuum to remove volatiles.



Scheme 2

To 4.00 mg of beads was added 2.0 ml of pyridine (M.W. = 79.10) and 100  $\mu\text{l}$  of trichloroethyl chloroformate (M.W. = 211.86). The mixture was stored at room temperature ( $25^\circ\text{C}$ ) for 24 hours. The mixture was filtered, washing five times with  $\text{CH}_2\text{Cl}_2$ . The beads and filter paper were put under vacuum to remove volatiles.

Individual beads were removed from each of the three sets of beads and placed on a Tacky Dot Array glass slide (from SPI Supplies). These glass slides have 40  $\mu\text{m}$  diameter adhesive dots, spaced 500  $\mu\text{m}$  apart. Individual beads were held on to the glass slide by placing them on the adhesive dots.

All samples were run on a Nicolet/Spectra-Tech microscope with an ATR objective with a ZnSe internal reflection element. This ATR objective can be used on any Nicolet/Spectra-Tech microscope. In addition, the InSpecIR™, which has a built-in ATR element, can also be used. All spectra were acquired with 8  $\text{cm}^{-1}$  spectral resolution with Happ-Genzel apodization, and 128 scans were co-added for each spectrum. The background spectra were obtained through air. The sample size was circular, with a diameter of 42  $\mu\text{m}$ .




Figure 2: Infrared ATR spectra of top) unreacted bead, middle) Scheme 1 reacted bead, and bottom) Scheme 2 reacted bead


## Results and Discussion

**Scheme 1:** The reaction product can be characterized from the infrared spectrum (middle spectrum, Figure 2). Spectral bands from the bead substrate can be seen (top spectrum, Figure 2), but new features have appeared in the spectrum of the reacted product. The appearance of a band at 1668  $\text{cm}^{-1}$  is indicative of the amide I stretch. The band at 1536  $\text{cm}^{-1}$  represents the amide II vibration. In addition, the growth of the bands at 1463  $\text{cm}^{-1}$  and at 1357  $\text{cm}^{-1}$  show an increase in the amount of  $\text{CH}_3$  present.


**Scheme 2:** The reaction product spectrum is the bottom spectrum shown in Figure 2. Again, spectral bands from the bead substrate are visible, but there are new bands present. The carbonyl stretch apparent at 1738  $\text{cm}^{-1}$  is indicative of the carbamate product. In addition, the appearance of a band at 1536  $\text{cm}^{-1}$  is due to the CNH group.

Conclusions - Continued on Page 24

# ElectroScan's Field Emission ESEM®




Magnetic Tape, 150,000x, 30kV, 5 Torr




Hole in resist, 60000x, 5kV, 3.5 Torr

# FE3030



Penicillium notatum mold, 8000x,  
5kV, 5 Torr



ElectroScan Corporation's ESEM FE3030, an ultra high resolution field emission gun Environmental Scanning Electron Microscope, designed, developed and manufactured jointly with Philips Electron Optics.

ElectroScan, ESEM, and Seeing Things You've Never Seen Before are registered trademarks of ElectroScan.

## Environmental Scanning Electron Microscope

- ◆ 2nm secondary electron resolution at 7 Torr at 30kV
- ◆ High stability Schottky Field Emission gun
- ◆ High beam current capability for x-ray on uncoated samples
- ◆ High resolution imaging of wet samples with no dehydration
- ◆ 45° Conical objective lens
- ◆ Motorized eucentric stage
- ◆ Microscope control from Windows®




**ElectroScan® Corporation**

Seeing Things You've Never Seen Before®

66 Concord Street  
Wilmington, MA 01887  
Tel (508) 988-0055  
Fax (508) 988-0062

Visit us on the World Wide Web—  
[http://www.tiac.net/users/  
scan/index.html](http://www.tiac.net/users/scan/index.html)

## NEW PRODUCT NEWS



### New LEICA GZ4 Stereomicroscope

The new LEICA GZ4 Stereomicroscope continues the Leica tradition of innovation at an economical price. In today's world of budgetary constraints, Leica provides the fundamental solution of value and quality stereoviewing with the LEICA GZ4.

Leica's Fundamental Solution:

- 4:1 zoom in a basic stereomicroscope with pre-centered optical mounting system assuring consistent alignment;
- 35 mm to 243 mm working distance for additional working space with specific applications;
- Magnification ranges from 2.1X to 180X with supplementary lenses;
- Designed with the unique, lightweight ZeroStat protective material which is specifically designed to dissipate static charge, preventing electrostatic discharge (EDS) and reducing dust attraction;
- A variety of stands is available to accommodate the

viewing of large specimens;

- Large zoom knobs provide a comfortable grip for even torque throughout the zoom motion.

Designed as a Greenough Optical System, the LEICA GZ4 provides high resolution and reduced eye fatigue without distortion or aberration over the entire field of view.

The LEICA GZ4 is manufactured at Leica's Optical Products Division, an ISO 9001 certified company. Leica strives to exceed customer expectations through individual commitment, teamwork and innovation. The LEICA is economically priced to meet challenging budget demands while maintaining high quality.

For further information regarding the LEICA GZ4 Stereomicroscope contact Leica, Inc., 111 Deer Lake Road, Deerfield, Illinois, 60015, Telephone: 800/248-0123, Fax: 847/405-0030

Circle Reader Inquiry #31

### Image Retriever Database Program Simplifies Multi-Media Image Management for Windows Users

Image Retriever, a new multi-media database for Windows, simplifies the capture, storage, retrieval, and printing of high-resolution image files with its intuitive user-interface and flexible database architecture. The innovative database program from Advanced DataBase Systems offers one-button image capture, OLE 2.0 drag-and-drop simplicity, advanced printing options, full-screen image slide show capabilities, and ODBC-compliance so users can query the database from other Windows programs.

Advanced DataBase Systems  
7931 South Broadway #322  
Littleton, CO 80122  
Tel.: (303)761-5635  
Fax: (303)761-3780

Circle Reader Inquiry #32

### 256 Gray Scale Printer has Windows Driver

Alden Electronics has announced the availability of a Windows driver for its Continuous Tone Printer, the 9315CTP.

The new software driver allows any Windows-based program, such as MS Powerpoint or Adobe Photoshop to output to the printer.

The 9315CTP provides 256 gray scale printing on a 10.1 inch-wide image for photographic quality hard copy. This quality makes the 9315CTP ideal for applications where a high quality image is needed. The printer uses economical thermal technology which keeps per-print costs as low as nine cents per piece. With just three moving parts, the printer has a Mean Time Between Failure of 20,000 hours, equivalent to about 10 years of full-time printing. It takes output from PC's, Macs, Sun and other workstations, and various imaging systems.

Alden Electronics  
40 Washington Street Tel.: (508)366-8851  
Westborough, MA 01581 Fax: (508)898-2427

Circle Reader Inquiry #33

### Direct Monitoring of Combinatorial Chemistry Reactions by Infrared Microspectroscopy - Continued from Page 22

#### Conclusion

Infrared microspectroscopy was used to detect and characterize reaction products present on resin bead substrates. ATR microspectroscopy was chosen because it requires only a single bead for analysis, and the technique is non-destructive to the sample. A bead is simply placed on a glass slide, contact is made between the bead and the ATR crystal, and a spectrum is obtained. Clearly, infrared microspectroscopy can enable combinatorial chemists to monitor their reactions at each step of the synthesis. Because there is no sample preparation involved, the entire process increases the overall efficiency of this new and exciting technique for organic synthesis.

1. "Special Report: Combinatorial Chemistry", Chemical & Engineering News, February 12, 1996, pages 28-73.
2. "Combinatorial Chemistry Speeds New Product Discovery by Michael Sullivan, Today's Chemist at Work, November 1995, pages 14-18.
3. "Infrared Spectrum of a Single Resin Bead for Real-Time Monitoring of Solid-Phase Reaction", by B. Yan *et al*, Journal of Organic Chemistry, Vol 60, 1995, pages 5736-5738

InspectIR is a trademark of Spectra-Tech, Inc.

### NORAN/Tracor Northern EDS Repair

Factory and field trained personnel with over 35 years total experience, located in the the Midwest, Southeast and on the West Coast, specializing in TN2000 and TN5500 repair, detector upgrades, data storage, imaging hardware and peripheral output devices - and used equipment resale.

Choose from a full maintenance or parts only contract, or on-demand service - at rates normally a fraction of others.

For further information, contact Doug Connors at:

### TN Analyzer Service

7897 Highway 19  
Dane, WI 53529

Tel.: (608)798-2005  
Fax: (608)798-1675