

DAB on Human Hepatocytes

Charles S. Gilbert, Carolinas Medical Center

Neutrophil myeloperoxidase standard 3,3'-Diaminobenzidine Tetrahydrochloride (DAB) procedure:

When staining neutrophils, a cell suspension would be preferred, but a finely minced buffy coat can be used, if you thick section and locate the stained cells before thin sectioning. Staining should be carried out as soon as possible after fixation. The myeloperoxidase (MPO) activity falls off such that we process samples within 2-3 weeks of fixation.

Note that fixation is for 1 hour. Enzyme activity may be destroyed if fixed for longer periods.

Chemical inhibition of some peroxidatic enzymes is possible, but inhibition is not 100% specific or effective. Chemical inhibition works best when an inhibition step (inhibitor and buffer) occurs before staining with the DAB solution which also contains the inhibiting agent.

The literature shows that it is difficult to assess different peroxidatic enzymes when they are colocalized using cytochemical techniques.

Always be sure to include appropriate controls (*i.e.*, a morphological control as well as cytochemical ones).

The DAB-stained cells are viewed without uranyl acetate-lead citrate (UALC) counterstaining, or with lead citrate (LC) staining only. Counterstains may confound the granule analysis.

We place our DAB incubation in the dark (in a cabinet).

This procedure stains neutrophil primary granules for myeloperoxidase activity [0.01% H₂O₂].

Procedure:

1) 3%. Glutaraldehyde in 0.1 M Cacodylate Buffer, pH 7.35 @ 40°C (on ice)

- 1 hour.
- 2) Rinse x 3 in 0.1 M Cacodylate Buffer + 7% Sucrose, pH 7.35. Store overnight in rinse buffer.
- 3) Rinse x 1 in DAB solution without H₂O₂.
- 4) DAB solution with H₂O₂ at room temperature 30 minutes. [Control: Use DAB solution without H₂O₂ for 30 minutes.]
- 5) Rinse x 1 in 0.05 M Tris Buffer, pH 7.6.
- 6) Rinse x 2 in 0.1 M Cacodylate Buffer + 7% Sucrose, pH 7.35.
- 7) 1% OsO₄ in 0.1 M Cacodylate Buffer, pH 7.35, at room temperature 1 hour.
- 8) Dehydrate and embed

Reagents:

DAB solution:

3,3'-Diaminobenzidine tetrahydrochloride (4 HCl) 5 mg.

0.05 M Tris Buffer pH 7.6 10 mL.

3% H₂O₂ 33 µL.

pH should be in the range of 7.0-7.6. When the buffer is fresh it is around 7.3-7.4.

DAB solution with 10 mM KCN:

3,3'-Diaminobenzidine 4 HCl 5 mg.

0.05 M Tris Buffer, pH 7.6. with 10 mM KCN 10 mL.

3% H₂O₂ 33 µL.

0.05 M Tris buffer, pH 7.6:

Tris-(hydroxymethyl)-aminomethane [99.0-99.5%] 6.06 gm.

make up to 1000 mL with distilled H₂O.

Adjust pH to 7.6.

3% H₂O₂ (Working Solution):

50% H₂O₂ (Commercial Stock Solution) X mL.

0.05 M Tris Buffer, pH 7.6 make up to 10 mL.

STOP HASSLING WITH MULTIPLE SERVICE CONTRACTS!

START by putting all of your instruments under one contract with MAS (regardless of make or model). Our expert EM SERVICE GROUP has the knowledge and skill to keep your instrumentation working at its best.

TEM'S/ SEM'S

**HITACHI
JEOL
AMRAY
CAMBRIDGE**

**TOPCON
ISI
ZEISS
PHILIPS**

We service all brands under one contract.

PREP EQUIPMENT

**ULTRAMICROTOMES
DUPONT & RMC
VACUUM COATERS
SPUTTER COATERS
MECHANICAL PUMPS
TURBO PUMPS
PLASMA ASHERS**

WE SERVICE

**COMPUTER CONTROLLERS
BEAM BLANKERS
WATER CHILLERS**

Contracts and On-Demand
Emergency Service at
Reasonable Rates from
Factory Trained Specialists.

1-800-421-8451

 **MAS**
MATERIALS ANALYTICAL SERVICES

3597 Parkway Lane • Suite 250 • Norcross, Georgia 30092 • 770-448-3200 • FAX 770-368-8256 or
616 Hutton Street • Suite 101 • Raleigh, North Carolina 27606 • 919-829-7041 • FAX 919-829-5518

ADVANCED ANALYTICAL PRODUCTS AND SERVICES

See "Verification" section below to determine the value of "X".

Note: 3% H₂O₂ working solutions should be made fresh before use using a concentrated stock solution which has been checked for the actual concentration. See the verification procedure for doing this spectrophotometrically.

10 mM potassium cyanide (KCN) [in 0.05 M Tris Buffer, pH 7.6]:

Potassium cyanide 65 mg.

Make up to 100 mL with 0.05 M Tris Buffer, pH 7.6

pH to 7.6 with 1 N HCl

MODIFICATIONS:

DAB with cyanide inhibition (10 mM KCN):

This procedure inhibits cytochrome oxidase and myeloperoxidase. 4-aminobenzhydrazide inhibits MPO. 1,2,3-Aminotriazole inhibits catalase.

For catalase localization use Tris buffer at pH 9.5-10.0.

For inhibition studies, always recheck the pH of the buffer-inhibitor solutions and pH them properly.

References for peroxidase (DAB) staining:

Angermuller, S. and H.D. Fahimi. 1981. Selective cytochemical localization of peroxidase, cytochrome oxidase and catalase in rat liver with 3,3'-diaminobenzidine. *Histochemistry* 71:33-44.

Biggar, W.D. and J.M. Sturgess. 1976. Peroxidase activity of alveolar macrophages. *Lab. Invest.* 34:31-42.

Breton-Gorius, J., Y. Coquin, and J. Guichard. 1978a. Cytochemical distinction between azurophils and catalase-containing granules in leukocytes. I. Studies in developing neutrophils and monocytes from patients with myeloperoxidase deficiency: Comparison with peroxidase-deficient chicken heterophils. *Lab. Invest.* 38:21-31.

Breton-Gorius, J. and J. Guichard. 1978b. Cytochemical distinction between azurophils and catalase-containing granules in leukocytes: Distribution in human promyelocytes and promonocytes. *J. Reticuloendothelial Soc.* 21:637-646.

Cotran, R.S. and M. Litt. 1970. Ultrastructural localization of horseradish peroxidase and endogenous peroxidase activity in guinea pig peritoneal macrophages. *J. Immunol.* 105:1536-1546.

Dunn, W.R., J. H. Hardin, S. S. Spicer. 1968. Ultrastructural localization of myeloperoxidase in human neutrophil and rabbit heterophil and eosinophil leukocytes. *Blood* 32:935-944.

Fahimi, H.D. 1969. Cytochemical localization of peroxidatic activity of catalase in rat hepatic microbodies (peroxisomes). *J. Cell. Biol.* 43:275-288.

Fahimi, H.D. 1979. An assessment of the DAB methods for cytochemical detection of catalase and peroxidase. *J. Histochem. Cytochem.* 27:1365-1366.

Graham, R.C., Jr. and M.J. Karnovsky. 1966. The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. *J. Histochem. Cytochem.* 14:291-302.

Kettle A. J., C.A. Gedy, M.B. Hampton, C.C. Winterbourn. 1995. Inhibition of myeloperoxidase by benzoic acid hydrazides. *Biochem. J.* 308:559-563.

Nishimura, E.T., G.M. Whest and H.Y. Yang. 1976. Ultrastructural localization of peroxidatic catalase in human peripheral blood leukocytes. *Lab. Invest.* 34:60-68.

Novikoff, A.B. and S. Goldfischer. 1968. Visualization of microbodies for light and electron microscopy. *Proc. Histochem. Soc.* 16:507 (Abstract).

Novikoff, A.B. and S. Goldfischer. 1969. Visualization of peroxisomes (microbodies) and mitochondria with diaminobenzidine. *J. Histochem. Cytochem.* 17(10):675-680.

Novikoff, A.B., P.M. Novikoff, N. Quintana, and D. Cleveland. 1972. Diffusion artifacts in 3,3'-diaminobenzidine cytochemistry. *J. Histochem. Cytochem.* 20:745-749.

Roels, F., E. Wisse, B. De Prest. and J. van der Meulen. 1975. Cytochemi-

cal discrimination between catalases and peroxidases using diaminobenzidine. *Histochemistry* 41:281-312.

Zellmer, D.M. and W.A. Shannon Jr. 1983. Morphometric cytochemistry of catalase and myeloperoxidase-containing granules in the rabbit polymorphonuclear leukocyte. *Histochem. J.* 15:211-230.

Verification of [H₂O₂]:

$$E_{240} = 43.6 \text{ m}^{-1} \text{ cm}^{-1} \text{ (extinction coefficient at wavelength 240)}$$

$$\text{MW} = 34.01 \text{ g/mole}$$

$$50\% \text{ H}_2\text{O}_2 = 1.2 \text{ g/mL (density, rho)}$$

Determine whether concentration is by weight or volume because serious miscalculations are possible.

Calculation for the molar concentration of:

$$50\% \text{ by wt} \Rightarrow (60\text{g})/120 \text{ g} = (60\text{g})/100 \text{ mL}$$

$$(60\text{g})/100 \text{ mL} \times 1/(34.01 \text{ g/mole}) = (1.761 \text{ mole})/100 \text{ mL} \Rightarrow 17.61 \text{ M}$$

Stock solution must be diluted to measure the Absorbance (A_{240})

l = length of light path in cm

C = concentration (moles/liter)

$$A_{240} = E_{240} l C \Rightarrow C = [A_{240} / (E_{240} l)] \times \text{dilution factor}$$

e.g., dilution of 1/1000 $A_{240} = 0.6375 \Rightarrow 14.6 \text{ M}$

$$(14.6 \text{ M} \times \text{MW})/10 = (49.6 \text{ g})/100 \text{ mL} = (49.6 \text{ g})/120 \text{ g} \Rightarrow 41.4\%$$

Use this percentage in making up the 3% H₂O₂ working solution.


Since the 50% H₂O₂ (17.61 M) was actually 41.4% (14.6 M) it follows that:

$$V_1 C_1 = V_2 C_2 \Rightarrow 14.6 \times = 1.06 (10 \text{ mL}) \Rightarrow x \text{ mL} = 10.6/14.6 = 0.728 \text{ mL}$$

Use 0.728 mL in making the 3% H₂O₂ working solution. ■

Light Element Peaks Revealed!

Oil Film on EDX Windows Removed:



Oily Window

Clean Window

Oil build-up on EDX detector windows can ruin sensitivity for light element X-rays in SEMs. To stop oil condensation and keep the system clean, smart SEM users rely on the XEI Scientific SEM-CLEAN™ system.

Result: Consistent light element X-ray results and contamination-free pictures. The Nitrogen purge of the inexpensive SEM-CLEAN system actively cleans your electron microscope while you're away.

SEM-CLEAN™ Stops the Oil

XEI
SCIENTIFIC

3124 Wessex Way, Redwood City, CA 94061-1348
650-369-0133 ■ Fax 650-363-1659

<http://www.msa.microscopy.com/SM/XEI/XEIHomePage.html>