


Powder Diffraction

K. E. Wiedemann, J. Unnam and R. K. Clark	Deconvolution of Powder Diffraction Spectra	130
K. E. Wiedemann, J. Unnam and R. K. Clark	Computer Program for Deconvoluting Powder Diffraction Spectra	137
R. J. Hill and I. C. Madsen	Data Collection Strategies for Constant Wavelength Rietveld Analysis	146
R. N. Rose	Charles S. Barrett: X-Ray Metallurgist — An Informal Portrait	163
G. Wilson and F. P. Glasser	Powder X-Ray Data for Synthetic Anhydrous Sodium Zirconium Silicates	176
Liu Dengfa and Wang Yuming	The "Hook Effect" of X-Ray Diffraction Peak Broadening of Multilayer Thin Films	180
G. I. Straguzzi and P. F. Dismore	Use of X-Ray Diffraction as an Analytical Technique for the Synthesis of Inorganic Catalysts	183
K. R. Waerstad and A. W. Frazier	X-Ray Powder Diffraction and Single-Crystal Data for the Isomorphous Series $(\text{Fe,Al})_3(\text{K,NH}_4,\text{H}_3\text{O})\text{H}_{14}(\text{PO}_4)_8 \cdot 4\text{H}_2\text{O}$	187
W. Wong-Ng <i>et al</i>	Standard X-Ray Diffraction Powder Patterns of Sixteen Ceramic Phases	191
	Departments	203
	Calendar of Meetings	204
	Meeting Reports	205
	JCPDS-ICDD Subcommittee Activities	206
	Short Courses and Workshops	207
	Book Reviews	208
	General Announcements	209
	Computer Comments	210
	New Products	212



Volume 2 Number 3 September 1987

Powder Diffraction An international journal of materials characterization

SIEMENS

Unequaled in x-ray diffraction

015.928

If one word can be used to describe x-ray diffraction as an analytical method, it has to be flexibility. No other technique demands the analysis of both the chemical and physical properties of materials, and no other method encompasses the broad range of sample forms.

The one word that describes the Siemens D 500 diffraction system is unequaled because it takes a diffraction system with unequaled flexibility to produce the results in today's laboratory environment.

The D 500 means high-speed diffractometry, 3-D graphics, automated phase analysis, texture or residual stress analysis, and gives you the proven performance and reliability to handle large volumes of samples regardless of their form or application.


D 500 — your unequaled choice.

Siemens AG, Analytical Systems E 689, D 7500 Karlsruhe 21, P.O. Box 21 1262, Federal Republic of Germany. Tel. (0721) 595-4295.

Siemens Energy & Automation, Inc., Analytical Systems, One Computer Drive, P.O. Box 5477, Cherry Hill, NJ 08034. Tel. (609) 424-9210.

Siemens . . . The company you can believe in

Set 37 Now Available
with Optional Level III Database



μ PDSM puts the world of XRD analysis at your finger tips.

Fein-Marquart has taken search/match out of the lab and literally put it at your finger tips with μ PDSM—a full featured software program that gives you mainframe power for your IBM PC. With both stand-alone and directly interfaced operation, μ PDSM provides cost-effective analysis on your PC. It delivers more accurate identification faster than any known alternative and gives you the 50,000+ pattern JCPDS file, minerals or common phases

subfiles, or your own library. Whether your x-ray diffractometer is manual, semiautomatic, or fully computer controlled, the μ PDSM is the answer to your direct interaction in solving simple or complex diffraction patterns. For more information or a presentation on μ PDSM, call or write FMA.





ERNEST F. FULLAM, INC.

Accessories for Microscopy • Scientific Consultants

900 Albany-Shaker Road • Latham, New York • 12110

Telephone (518) 785-5533

DIFFRACTION PATTERN MEASUREMENT



No. 19301: EFFA DIFFRACTION PATTERN MEASURING DEVICE

- Allows rapid and accurate measurement of "d" spacings from electron and X-ray diffraction patterns.
- Produces a chart record of variations in film density along a selected track.
- Line spacing can be read to 0.001mm for "d" spacing determinations.
- Extremely weak lines, line segments and spots may be measured.
- Handles all standard film and plate sizes on which diffraction patterns are recorded.
- Three different scan rates.
- Accurate measurement of angular relationships in spot patterns.
- Bench size required: 44" x 16".

Write or call for full catalog of Accessories for Microscopy

CANADA:
Micro Biological Supplies
41 Maple Avenue
Richmond Hill, Ont. L4C 6P4, Canada

CONTINENTAL EUROPE:
Touzart & Matignon
8 Rue Eugene Henaff
94403 Vitry-sur-Seine, France

GREAT BRITAIN:
Graticules, Ltd.
Morley Road
Tonbridge, Kent TN9 1RN England

Powder Diffraction

An International Journal of Materials Characterization

Deane K. Smith

Editor in Chief
Dept. of Geosciences
The Pennsylvania State University
239 Deike Building
University Park
Pennsylvania 16802 USA

J. L. de Vries

European Editor
Humperdinklaan 51
Eindhoven 5654 PB, Netherlands

Richard N. Rose

Assistant Editor and
Manager of Publication

Editorial Advisory Board

A. Albinati, Milan, Italy
L.V. Azaroff, Storrs, Connecticut
C. Baerlocher, Zurich, Switzerland
C.S. Barrett, Denver, Colorado
P. Bayliss, Calgary, Alta., Canada
C.Z. Bojarski, Katowice, Poland
A. Brown, Nykoping, Sweden
L.D. Calvert, Melbourne, Australia
D. Cox, Upton, New York
W. Eysel, Heidelberg, West Germany
J. Fiala, Pizen, Czechoslovakia
V.A. Frank-Kamenetsky, Leningrad, U.S.S.R.
L. Frevel, Midland, Michigan
P. Gado, Budapest, Hungary
H. Goebel, Munchen, West Germany
Th. Hahn, Aachen, West Germany
G.G. Johnson, Jr., State College, Pennsylvania
Q. Johnson, Livermore, California
J.I. Langford, Birmingham, U.K.

Ron Jenkins

Managing Editor
JCPDS-International Centre
for Diffraction Data
1601 Park Lane
Swarthmore
Pennsylvania 19081 USA

Mary M. Rossi

Assistant to the
Managing Editor

Ron Anderson

Departments Editor
Powder Diffraction
IBM Corporation
Dept. 13W, Bldg. 630-E70
Hopewell Junction
New York 12533

Gregory J. McCarthy

Editor for New Diffraction Data
North Dakota State University
Department of Chemistry
 Fargo, North Dakota 58105-5516 USA

Eloise Humez Evans

Editorial Assistant

D. Louër, Rennes, France
H.F. McMurdie, Washington, District of Columbia
M.E. Mrose, Washington, District of Columbia
M.H. Mueller, Argonne, Illinois
M. Nichols, Livermore, California
B.H. O'Connor, Bentley, Australia
W. Parrish, San Jose, California
B. Post, Brooklyn, New York
E. Prince, I.U.Cr. Representative
R. Shirley, Guildford, U.K.
R. L. Snyder, Alfred, New York
W. Su, Changchun, China
P. Suortti, Helsinki, Finland
J.W. Visser, Delft, Netherlands
S. Weissmann, Piscataway, New Jersey
A.J.C. Wilson, Cambridge, U.K.
T. Yamanaka, Tokyo, Japan
R.A. Young, Atlanta, Georgia
L. Zevin, Beer-Sheva, Israel

Publisher

JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081, U.S.A.

Powder Diffraction is a journal of practical technique, publishing articles relating to the widest range of application — from mineral analysis to epitaxial growth of thin films and to the latest advances in software. Although practice will be emphasized, theory will not be neglected, especially as its discussion will relate to better understanding of technique.

Powder Diffraction is published four times annually by the JCPDS-International Centre for Diffraction Data.

Manuscript submissions. The Editors will consider all manuscripts received, but assume no responsibility regarding them. Materials will be returned only when accompanied by appropriate postage.

Subscriptions. The annual subscription rate in the United States and Canada is \$39.50; Library: \$75.00; other than U.S.A., Canada and the Far East, the annual subscription is \$55.00. Subscriptions to the Far East, including Japan, China, Taiwan, Malaysia, the Philippines, Indonesia and Korea should be made via Sanyo Shuppan, Boeki Co., Inc., P.O.B. 5037, Tokyo International 100-31, Japan. Airmail delivery available for subscribers outside U.S./Canada for an additional cost of \$35.00 per volume (4 issues).


Payment may be made in U.S. dollars by company check, or international money order. Please address communications to the publisher's office.

Advertising. For advertising rates and schedules contact the Publisher's Office JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, PA 19081, Telephone (215) 328-9405.

Reprints and permissions. Contact the Publisher's office.


Postal Information. Powder Diffraction (ISSN 0885-7156) is published quarterly for \$39.50 a year (U.S. and Canada) by the JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081. JCPDS principal office: 1601 Park Lane, Swarthmore, Pa. 19081. Julian Messick, Jr., General Manager. © 1987 JCPDS-International Centre for Diffraction Data. Postmaster: Send address changes to JCPDS-International Centre for Diffraction Data, 1601 Park Lane, Swarthmore, Pennsylvania 19081.

International CODEN Service (Intercode): PODIEZ
ISSN 0885-7156
Telex: 847170



Real Simultaneous X-Ray Diffraction Is Here.

Data Collection Increased Up To 500%
Over Conventional Scanning Techniques.



The Inel CPS Diffraction System.


In Powder Analysis:

- No movement of detector or the sample.
- High-speed Data Collection displayed in Real Time.
- Extremely versatile horizontal or vertical configuration, addition of high temperature and other accessories.
- The Detector is easily adaptable to most existing Diffraction Systems.

In other applications:

- Texture analysis with simultaneous generation of 3 pole figures.
- Can be used for Single Crystal, Stress, Amorphous and Kinetic studies.

3-Dimensional Graphic Output of a Kinetic Study.



American Instruments, Inc. 201-636-5770

Telex: 5101006499 AII NJ 185 Port Reading Avenue, Port Reading, N.J. 07064

American Instruments features State-Of-The-Art Inel Products.

inel - 1261, rue L. Blériot, 78530 Buc - France Tél. (33)-1 39.56.31.90 - Télex 698502 F

Powder diffraction is our specialty.

In fact, it's all we do.


That's why we can deliver all those technological advances others have been promising. And we give you the largest selection of software there is. That's also why Scintag users enjoy the kind of after-sale support that helps them realize substantial returns on their investment, year after year.

Take our solid-state detector as an example. It's routinely available on the PAD V system, and it works for both horizontal and vertical configurations. We were the first to offer time averaged and signal processing data collection schemes for routine analysis. No one else comes close to our true microstep data collection technology ($1^\circ=3200$ steps). Our continuously variable radius goniometer remains unique.

Peruse our software offerings, and you'll see how we've translated our powder diffraction expertise into practical solutions. In packages that solve everyday problems and address specialized needs for metallographic, crystallographic and thin film applications.

To help our users apply this advanced technology to their analytical problems, Scintag maintains the best customer support program in the industry. Through a combination of applications assistance and fast service turnaround, our users are able to fully exploit the capabilities of our systems.

Contact the powder diffraction specialists today. Better yet, ask for a demonstration, and see for yourself.



Scintag, Inc.

3350 Scott Boulevard, Building 28
Santa Clara, CA 95054
(408) 748-8544 TLX 171132

PD5

Finally. An affordable, easy-to-operate alternative to high cost diffractometer automation.

The Databox.

The acquisition package:

The Databox is a stepping motor driver and data collection system specifically designed to control your diffractometer, all in the space of a two-wide NIM module. Just talk to it from a computer terminal or your PC using an incredibly friendly command language, and all your data acquisition needs are satisfied.

The analysis package:



You can also buy the Databox bundled with MDI's Peak Identification and Micro-ID Search/Match software to run on your IBM PC (or compatible), giving you a complete x-ray control, acquisition, and analysis system.

The bottom line:

For well under \$5000, the Databox will fully automate your x-ray equipment. Add the analysis software, and the total cost is just over \$10,000.



Finally. A system with proven reliability and performance at a reasonable price. To learn more, contact us at:

Radix Instruments
1019 Stratford Ave.
South Pasadena, CA 91030
(818) 441-5351




Materials Data, Inc. has more XRD Software for your IBM-PC

Micro-Peak performs data reduction including smoothing, background subtraction, peak finding, graphical comparisons and more. Combine with **Micro-ID** for a complete pattern analysis system.




Micro-Powd is the latest and best pattern simulation program. Based on 20 years of mainframe use. Now available for your PC.

Micro-ID provides fast, effective Search/Match to identify unknown materials from powder patterns. Use JCPDS files (full file, subfiles, CD-ROM) or your own data.



Micro-Laue




Micro-Shadow enables you to do profile fitting with a variety of functions. This pattern analysis system includes features for size and strain analysis as well as diffractometer calibration.

Materials Data, Inc. specializes in PC software for the X-ray laboratory. Other programs are available. For example, **Micro-Laue** is a program to simulate front- or back-reflection Laue patterns. Contact us concerning your software needs. Our telephone number is (415) 449-1084.

MATERIALS DATA, INC.
P.O. Box 791
Livermore, CA 94550
(415) 449-1084

SIEMENS

The Choice You Can Believe In The D 500 X-ray Diffraction System



X-ray diffraction has become one of the most important and necessary tools for solving analytical and research problems in laboratories throughout the world, and the Siemens D 500 diffractometer provides the universal capabilities for satisfying the widest range of applications.

The D 500 is the one system that combines a wide selection of accessories with a comprehensive but user-oriented software package that enables you to perform routine powder diffraction analysis in production, control or sophisticated research by the addition of options, such as stress, texture, high-temperature, and Theta-Theta diffractometry — or dramatically speed up the analysis with

the addition of a position sensitive detector.

Already operating in numerous laboratories worldwide, the D 500 is evidence of Siemens commitment to satisfying all your present and future needs in x-ray diffraction. Siemens is the choice that you can believe in for x-ray diffraction systems. For more information write or call: Siemens AG, Analytical Systems E 689, D 7500 Karlsruhe 21, P.O. Box 21 1262, Federal Republic of Germany, Tel: (0721) 595-2425, Telex: 78255-69. Siemens-Allis, Inc., Analytical Systems, One Computer Drive, P.O. Box 5477, Cherry Hill, NJ 08034, (609) 424-9210.

Siemens . . . The company you can believe in.


PD8

PC-PDF


The Powder
Diffraction File
on CD-ROM



Search on key
fields within
the data base.



Search on
combinations
of fields.




Display results within seconds.

PC-PDF, a low cost personal computer based system, makes the entire Powder Diffraction File available on a single CD-ROM disk ... and, through use of optimum packing and access algorithms, displays results within seconds.

- Use with an IBM-XT or equivalent system.
- Annually updated disks will be made from the full data base and will include all additions and corrections.
- Space saving: One disk contains data equivalent to a file cabinet full of Data Cards.



You can now enhance your present system for search, retrieval and display.
For additional descriptive material and ordering information, contact:

JCPDS-International Centre for Diffraction Data
1601 Park Lane, Swarthmore, PA 19081-2389 USA
(215)328-9400 Telex: 847170



COMPONENTS AND SYSTEMS FOR SEMICONDUCTOR AND MATERIALS CHARACTERIZATION...

BLAKE INDUSTRIES



Blake Industries is the leading manufacturer of x-ray components and systems that meet the requirements of modern materials and semiconductor technology.

We develop—and constantly evolve—instrumentation sensitive to new materials and to the demands imposed by smaller and smaller semiconductor circuits. Our product designs reflect today's major emphasis on detecting material defects, as well as the need for analyzing epitaxial films by measuring lattice mismatch between the epitaxy and the film. And our instrumentation is known as the

most flexible in handling a wide range of materials with all the necessary radiation safety.

Other leading companies, such as IBM, Bell Labs, Westinghouse, TRW and AT&T, are putting Blake Industries' technology to work for them. Isn't it time yours did?

BLAKE INDUSTRIES, INC.


660 Jerusalem Road
Scotch Plains, NJ 07076
(201) 233-7240
Telex: 7109979581



PD10

A New Generation Of X-Ray Generators

Spellman, a leader in high voltage technology, has developed a power supply with the features diffraction scientists require. The DXR 3000 offers a unique combination of low cost, high performance, and compact size. These features are achieved through Spellman's experience in high frequency resonant inverter technology.



FEATURES:


- Greater than 90% efficiency.
- High stability through precision feedback control circuits.
- Solid encapsulants insure maintenance-free operation.
- System fault diagnostics and display, including arc detection, filament short, overvoltage, and tube water flow interlock.
- Automatic slew of the high voltage and emission current to preset values.
- Current source provides automatic matching of the filament current to X-Ray tubes with various heater voltages.
- Continuous adjustment for both high voltage and emission current.

TECHNICAL DATA: MODEL DXR3000

Output Power	3KW
High Voltage Output	0-60KV
Emission Current	0-50mA
Emission Current Regulation	± 0.1% for a 50% voltage change
Emission Current Accuracy	± 1%
Line Regulation	± 0.01%
Ripple	0.1% RMS
Input Voltage	220V ± 10%/50 or 60Hz single phase
Input Power	4.8KVA
Filament Output	0-6A/20VDC Compliance
Dimensions	Three chassis each 7"H × 19"W × 20"D
Weight	180 lbs total

OPTIONS:

- 6KW and 9KW available.
- Outputs up to 100KV available.
- 208V 3 phase model available.
- Remote interface for all analog and digital functions available.
- Fast emission current response (stable to 0.1% within 50mS) available.



Spellman HIGH VOLTAGE ELECTRONICS CORP.

Our Fourth Decade of Leadership in High Voltage Technology.

For more information please write or phone: Spellman High Voltage Electronics Corporation
7 Fairchild Avenue, Plainview, New York 11803 • (516) 349-8686 • TWX510-221-2155

PD11

X-ray Diffraction Reference Standards and Zero-background Sample Plates Custom Designed and Built for any Application

Your first step to improved x-ray diffraction results should be to contact The Gem Dugout for quality diffraction alignment standards and zero-background plates. And the next step is successful x-ray diffraction results.

The Gem Dugout
1652 Princeton Drive
State College, PA 16803
(814) 865-5782

PD12



DIANO X-RAY ANALYSIS MAXIMUM PERFORMANCE AT MINIMUM COST

Whatever your process control or analysis application – chemicals, metals or the environment itself, one source offers highest performance at lowest cost — DIANO. Upgradable design and skilled field service assure that more than 2,000 diffractometers and spectrometers worldwide continue a 30-year record of maximum performance.

A DIANO system is a long-term investment that offers immediate improvement in accuracy, downtime and productivity.

DEMAND DIANO – Technology That Never Goes Out Of Style.



DIANO 30 Commerce Way, Woburn, MA 01801 U.S.A.
Tel. (617) 935-4310 • Telex 94-9306

PD13

Requiem for a heavyweight

The fact is, our new Psi Peltier cooled silicon detector obsoletes monochromators and proportional detectors and their "blind" PHA electronics.

Psi is small and lightweight, mobile and worry free—the first to deliver solid state detector performance without the clumsiness of liquid nitrogen.


With Psi, you get 3 to 5 times faster throughput, 10 times better energy resolution. And the best diffraction patterns you've ever seen. Psi is like having an electronic monochromator. You get clear $K\alpha$ or $K\beta$ patterns from any tube anode.

Teamed with spectrum analyzer electronics, Psi gives you the complete sample spectrum. Fully resolved. This means you choose your peaks, then optimize the patterns to fit your application—Compton-Rayleigh for instance. Or XRF/XRD studies.

Better still, Psi comes with mechanical and electronic interfaces that make it 100-percent compatible with existing diffractometers. Not just new ones.

For a list of leading companies now using Psi with their systems—plus complete technical and upgrade information—write or call toll-free today. 1-800-227-0277.

Psi. A technical knockout.




Kevex Corporation 1101 Chess Drive, Foster City, CA 94404 415-573-5866 1-800-227-0277
Kevex Affiliate (Getac GmbH) Tel 49-6131-40091 Kevex Canada, Ltd. Tel 416-731-2161 Kevex France (Fondis, S.A.) Tel 33-1-34810024 Kevex UK Tel 44-582-400596



PD14

Rocklabs Laboratory Ring Mills

The Proven Pulveriser



Two Way Air Control Valve
Clamp Dust Cover
Clamp Mounting Plate
Flexible Air Line
Laminated Timber Noise
Reduction Cabinet



Rocklabs Laboratory Ring Mills have proven their superior pulverising capabilities since 1969 in 400 laboratories throughout the world. Consider Rocklab's unique quality features and you'll understand why!

- Largest range of heads in the world — one gram to 1000 grams in steel, tungsten carbide, ceramic and agate.
- All heads incorporate one-piece bowls with rounded corners for maximum pulverising efficiency, easy cleaning and long maintenance-free life.
- Pneumatic or mechanical clamping options.
- All heads are made of high purity materials, all exclusive to Rocklabs, for lowest sample contamination.
- Noise reduction enclosure, to better the international 85 dB maximum noise level, included in price.
- Competitive Pricing.
- Urgent replacement parts service.

ROCKLABS

For further information on Rocklabs Laboratory Ring Mills, contact Rocklabs PO Box 18-142 Auckland 6, New Zealand. Telephone (09) 574-698. Telegraphic "Rocklabs" Telex NZ 60550 "ROCKLAB". Agents in Australia, Canada, Chile, Indonesia, Peru, Philippines, UK, USA.

PD15

PORTABLE IMAGE X-RAY INTENSIFIER



FEATURING . . .

- REAL TIME CAPABILITIES FOR . . .
X-RAY IMAGING
X-RAY TOPOGRAPHY
LAUE CRYSTAL ORIENTATIONS
- RAPID ALIGNMENT OF . . .
CONVENTIONAL TOPO CAMERAS
COMPLEX X-RAY OPTICS


OFFERING . . .

- DIGITAL IMAGING
- INTERFACE/IBM PC
- CUSTOM DESIGN SERVICES



brimrose corporation of america
baltimore, maryland 21236
301/668-5800 • telex: 910-997-6817

PD16



Simultaneous detection using position-sensitive proportional counter (PSPC) . . .

- Data acquisition 100 times faster than classical step/scans
- Optimum detector resolution with exclusive low-noise, carbon-coated quartz anode
- No gas support required with sealed counter design
- Available fully integrated with NIM-based signal processing and PC-based multichannel analyzer

Researchers and engineers performing small angle x-ray scattering, time resolved diffraction, line shape or crystalline structure analysis will benefit from TEC's PSPC-based low energy x-ray diffraction system. For more information, call (615) 966-5856.




TECHNOLOGY for ENERGY CORPORATION
LEXINGTON DRIVE P.O. BOX 22996 KNOXVILLE, TN 37933-0966 USA

PD17

A new class of diffractometer with X-RED

X RED... an X Ray Electronic Diffractometer based on a Curved Position Sensitive detector ranging 120 degrees with a resolution of better than 0.03 dg. and its associated electronic, X-ray generator, water refrigerator unit, all included in an independent X-ray proof cabinet. Of course linked to a PC Olivetti Computer with its fine graphic performance, to display in real time a diffraction pattern. And more, automatic sample exchanger for capillaries, liquid, high, low temperature attachment. X-RED - CAN SAVE YOU A LOT OF TIME IN DATA COLLECTION SPECIALLY IN THE FIELD OF ROUTINE POWDER AND AMORPHOUS MATERIALS ANALYSIS.




38066 RIVA DEL GARDA - Zona Industriale Baltera (ITALY) - Tel. 0464/553426 - Telex 400278 RIVTOUR

PD18

Supper

Complete Debye-Scherrer Powder Camera Systems Accessories, and Film




114.6mm Powder Camera with optional viewing stand, illuminator, and Gandolfi "Single Crystal Randomizer" which is interchangeable with standard powder sample holders

Two Sizes: 114.6mm camera, 57.3 camera
Easy to load and unload in the darkroom, the cameras are precision manufactured with no screw-on components. The collimator and beamtrap are accurately secured in position magnetically and a removable sample holder permits insertion of the specimen sample outside the body of the camera.

Guinier Camera for transmission and reflection photographs is also available.

Simplified Specimen Line-up and Darkroom Handling Procedures Offer the Ultimate in Convenience

- Gandolfi "Single Crystal Randomizer"
- Sturdy Viewing Stand and Illuminator
- Precision Film Punch and Trimmer
- Universal Track and Tripod Mount
- Large Inventory of Capillary Tubes
- Large Inventory of Kodak 35mm Film
- Film Measuring Device and Illuminator



Reads accurately to 0.01mm and eliminates eyestrain

Manufacturers of Fine X-Ray Diffraction Equipment For Over 45 Years


Charles Supper Company
INCORPORATED

15 TECH CIRCLE, NATICK, MA 01760, U.S.A.
TELEPHONE: (617) 655-4610

PD19

THE WAY OF THE FUTURE FOR X RAY DIFFRACTION

AMORPHOUS




SINGLE
CRYSTAL

TEXTURE

KINETIC
STUDIES

STRESS

PHASES
ANALYSIS



OVER 120 DEGREES !

CPS 120 SYSTEM BRINGS THE FUTURE IN YOUR
LABORATORY

ULTRA FAST - ACCURATE - RELIABLE

THE CPS 120 SYSTEM - CURVED POSITION SENSITIVE DETECTOR

ASK FOR MORE INFORMATIONS

inel

1261, rue L.-Blériot, 78530 BUC - FRANCE - Tél. (1) 39 56 31 90 - Télex 698502 F

U.S.A.

AMERICAN INSTRUMENTS, Inc.
185 Port Reading Ave.
PORT READING
N.J. 07064
Téléphone : (201) 636-5770
Télex : 5101006499


ITALY

ITALSTRUCTURES
38066 Riva del Garda
Zona Industriale Baltera
Téléphone : 0464553426
Télex : 400278

GREAT BRITAIN

ANASPEC INTERNATIONAL LIMITED
Anaspec House, Faraday Road
NEWBURY
BERKSHIRE RG 13 2 AD
Téléphone : 0635 35733
Télex : 847509 (Cables optics G)

PD20




Nicolet X-ray: Integrated solutions to the diffraction puzzle.

Nicolet X-ray brings the full capability of X-ray diffraction to the polycrystalline laboratory. Whether you need a versatile powder diffractometer system or a specialized system for materials analysis, Nicolet's comprehensive line of powder diffractometers and advanced technology attachments can provide your laboratory the optimal blend of research performance and experimental versatility. Our instruments range from biplaner goniometers to transmission systems to double monochromator diffractometers. In addition, Nicolet systems offer leading technology attachments and accessories.

And now the uncompromising combination: Nicolet's new x-ray powder diffractometer with the MicroVAX 2000™ computer.

With Nicolet's new I2/V-2000, you no longer have to sacrifice powder diffractometer performance and versatility for computer compatibility. Incorporating the impressively compact MicroVAX 2000 32-bit architecture computer, the I2/V-2000 offers versatile application performance and unsurpassed processing speed and power for today's operating environments. The system can also be linked into the Digital Equipment Corporation LIMS (Laboratory Information Management System) network.

The I2/V-2000 is manufactured in the U.S.A. and supported internationally by more than 200 Nicolet Sales and Service representatives worldwide. Ongoing training and service programs are available. For a copy of our new product brochure or to schedule a demonstration, contact:
Nicolet Analytical Instruments 5225-1 Verona Road,
Madison, WI, 53711, Phone 608/271-3333.



digital™ Cooperative
Marketing
Program

MicroVAX 2000 is a registered trademark of
Digital Equipment Corporation

PD21 For Information
PD22 For Demonstration

Nicolet
INSTRUMENTS OF DISCOVERY

Editorial

X-Ray Wavelength Nomenclature

There has been increasing use of X-ray methods for materials characterization since the early 1950's, and the initial use of the diffraction and emission methods goes back four decades before that time. Along with this increase in the application of X-ray methods there has been a growing problem with the system of X-ray wavelength nomenclature. The current system was proposed by Siegbahn back in the 1920's and since his initial proposal, more lines have been identified, many of which do not fit with the nomenclature system. Over the years, this system has led to an increasing number of inconsistencies within a given line series. As an example, in the K series, the $3p \rightarrow 1s$ transitions are labelled $\beta_{1,3}$ rather than $\beta_{1,2}$, because it was originally believed that the line was a singlet, along with the $4p \rightarrow 1s$ transition also believed to be a singlet. Much later, as the resolution of crystal spectrometers improved, it was found that both of these lines were doublets. Because the $4p \rightarrow 1s$ transition was already called β_2 , and because β_3 was already allocated to the $3p \rightarrow 1s(1/2)$, the lines in the newly found doublet were called $K\beta'$ and $K\beta''$.

With the identification of further lines, especially in the longer wavelength series, and with the development of the Auger and other electron spectroscopic techniques, this problem has become intolerable. As an example, a reasonably modern X-ray spectrometer is well able to detect lines from the M series of all elements from the rare earths and on up in the periodic table. At this time, less than one half of these lines have been allocated names within the Siegbahn system!

Commission V.4 of the International Union of Pure and Applied Chemistry has for more than 20 years been addressing the problems of nomenclature, symbols and their usage in the various analytical spectroscopic disciplines.

This group has recently recommended¹ a brand new notation for X-ray spectroscopic wavelengths to be called "IUPAC NOTATION". The proposed system is based on energy level designations and has the advantage of being simple, logical and easy to apply to all types of transition. Moreover, it is consistent with the notations used in electron spectroscopy. As well as being easier to apply, the new system has the added advantage of relying solely on the ASCII character set. With the increasing use of automated typesetting machines for drawing up wavelength tables, there should be no further confusion caused by the inability of the typesetter to print certain characters. The X-ray levels are defined as those states which occur as initial or final states in X-ray transitions. The X-ray level notation follows earlier conventions except for the use of Arabic rather than Roman numerals. As the main rule, transitions between X-ray levels are denoted by the level symbols for the initial and final state separated by a hyphen. The initial state is placed first, irrespective of the energetic ordering. Thus, the $K\alpha_1$ and $K\alpha_2$ will now be known as the K-L3 and K-L2 respectively.

The Editors of *Powder Diffraction* are now seeking permission from IUPAC to reproduce the full report of the new wavelength nomenclature system in a future issue. It is hoped that the new system will soon be ratified, and it should prove a much more practical solution for X-ray line nomenclature.

1. Jenkins, R., Manne, R., Robin, J. and Senemaud, C., "Provisional Nomenclature Report from Commission V.4 — Nomenclature, Symbols, Units and their usage in spectrochemical analysis; Part VIII: Nomenclature system for X-ray spectroscopy", IUPAC (1987).

Ron Jenkins
Principal Scientist
International Centre for Diffraction Data