

THE JOURNAL OF HYGIENE

(Founded in 1901 by G. H. F. Nuttall, F.R.S.)

EDITED BY

R. M. FRY, M.A., M.R.C.S.

IN CONJUNCTION WITH

R. D. BARRY, M.A., B.V.Sc., Ph.D.

W. I. B. BEVERIDGE, D.V.Sc., M.A.

A. W. DOWNIE, M.D., F.R.S.

R. I. N. GREAVES, M.D., M.R.C.P.

J. O. IRWIN, M.A., D.Sc., Sc.D., F.S.S.

E. T. C. SPOONER, M.D., F.R.C.P.

Sir GRAHAM WILSON, M.D., LL.D., F.R.C.P., D.P.H.



CAMBRIDGE UNIVERSITY PRESS

Bentley House, 200 Euston Road, London, N.W.1


American Branch: 32 East 57th Street, New York, N.Y.10022

Subscription per volume £4. 10s. (\$17.50) (4 parts)

Single parts 32s. 6d. (\$5.00)

A NEW PURIFIED AGAR FOR IMMUNO-ELECTROPHORESIS

Photograph by courtesy of Department of Clinical Photography, St. James' Hospital, London



Immuno-Electrophoresis performed on Oxoid I.D. Agar using Oxoid barbitone acetate buffer BR 11g.

- 1 Serum from a case of myelomatosis. Note the deformed IgG, particularly the cathodal end which has reacted only poorly with the antiserum and shows a broad, boat-shaped widening.
- 2 Serum from a case of biliary cirrhosis (primary). The pattern resembles that from a normal serum except that all three immunoglobulins are increased and are readily visible. IgG is the only immunoglobulin seen in normal patterns. IgA and IgM, here well marked, are not normally seen at all.

Alb	Albumin	Cp	Ceruloplasmin
Om	Orosomucoid	Hp	Hemopexin
Ht	Haptoglobin	α_2M	α_2 -macroglobulin
Lp	β -lipoprotein	Tf	Transferrin
IgA	Immunoglobulin A (Syn: γ_1A ; β_2A)		
IgM	Immunoglobulin M (Syn: γ -macroglobulin)		
IgG	Immunoglobulin G (Syn: γ_2 ; 7S γ)		
IgG[M]	Pathological IgG (Myelomatosis)		
β_1A	} Subfractions of the lipoprotein		
β_1C			

Oxoid I.D. Agar has been developed especially to meet the exacting requirements of Immuno-Electrophoresis. This technique demands an Agar of low conductivity, virtually free from chloride and sulphate ions. All pre-treatment, washing, etc., is eliminated.

CONVENIENT TABLET FORM Oxoid I.D. Agar comes in convenient tablet form and each tablet makes 50 ml. of Agar gel.

For a leaflet on the technique of Immuno-Electrophoresis on Agar and a comprehensive manual containing details of the Oxoid range, write to: OXOID LIMITED, SOUTHWARK BRIDGE ROAD, LONDON, S.E.1.



for progressive laboratories