






Microscopy TODAY

Volume 30 Number 6 2022 November





Wood sample scanned at 280 nm voxel size - vessels are color-coded to thickness.

TESCAN UniTOM HR

The first micro-CT system to provide sub-micron spatial resolution and high temporal resolution dynamic CT in a single, highly versatile system.

**What can you do with TESCOAN UniTOM HR?
Contact us today to find out:**



Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in light, electron, or scanning probe microscopy, or the biological or the physical sciences, MSA takes your knowledge to the next level!

Members Receive:

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and more!



Join MSA Today!
Visit www.microscopy.org

DiATOME Diamond Knives

the incomparable
Diamond Knife for
all fields of research...

Please visit our website
for our complete range
of Diamond Knives:
www.diatomeknives.com

cryo AFM, ultra AFM

*Extremely smooth sample surfaces
Best possible structure preservation*

These high-quality DiATOME diamond knives are not just perfect for producing ultrathin and semi-thin sections, but also for surface sectioning of all kinds of biological and industrial specimens for AFM investigation.

Instead of a section as for TEM, the specimen surface is investigated using AFM.

Morphology of a blend of two SBS block copolymers with different chain-architecture. AFM tapping mode, phase image, image size = 3x3µm. Rameshwar Adhikari, Institut für Werkstoffwissenschaft, Martin-Luther-Universität, Halle-Wittenberg.

In order to achieve the best results for AFM investigation, use the highest quality diamond knives.

DiATOME ultra AFM and cryo AFM knives are specially tested to ensure that they meet the increased quality requirements of AFM investigation. They produce extremely smooth sample surfaces and guarantee the best possible structure preservation. Available in 35° with a cutting range of 10-100 nm in 2 mm and 3 mm sizes.




trimtool 20, 45, 90

One trimming tool for all of your trimming needs, be it at room or cryo temperatures.

A well trimmed sample is a precondition for perfect section ribbons. Due to the extreme sharpness of our diamond blades, less mechanical damage is applied to the sample during trimming. Very shiny sample faces and precise sides are the result.

AFM amplitude image of the muscle of cat's mite *Otodectes cynotis*. The contrast covers amplitude variation in the 1-3nm range. Size of the whole image equals 4.6microns. Nadejda Borisovna Matsko, Institut für angewandte Physik, ETH Zürich.




120,000x magnification. Ronald Walter, BASF Aktiengesellschaft, Polymer Physics, D-67056 Ludwigshafen.

DiATOME U.S.

P.O. Box 410 • 1560 Industry Rd. • Hatfield, Pa 19440
Tel: (215) 412-8390 • Fax: (215) 412-8450
email: info@diatomeknives.com • www.diatomeknives.com

Contents

About the Cover



Fourier ptychography process for analyzing images. For further information please see Aidukas et al. on pages 40–45.

TEM Development

- 10 A Roadmap for Edge Computing Enabled Automated Multidimensional Transmission Electron Microscopy**
 Debangshu Mukherjee, Kevin M. Roccapriore, Anees Al-Najjar, Ayana Ghosh, Jacob D. Hinkle, Andrew R. Lupini, Rama K. Vasudevan, Sergei V. Kalinin, Olga S. Ovchinnikova, Maxim A. Ziatdinov, and Nageswara S. Rao
- 20 An Interactive Simulation and Visualization Tool for Conventional and Aberration-Corrected Transmission Electron Microscopy**
 X.-Z. Li

X-Ray Microscopy

- 28 Achieving High-Resolution Hard X-ray Microscopy using Monolithic 2D Multilayer Laue Lenses**
 Wei Xu, Weihe Xu, Nathalie Bouet, Juan Zhou, Hanfei Yan, Xiaojing Huang, Ming Lu, Maxim Zalalutdinov, Yong S. Chu, and Evgeny Nazaretski

Biological Applications

- 34 Ultraviolet Autofluorescence Microscopy of *Nanotyrannus lancensis* Sections Reveals Blood Clots in Vessel Canals**
 Mark H. Armitage

Microscopy Education

- 40 Applications and Extensions of Fourier Ptychography**
 Tomas Aidukas, Felix Wechsler, Lars Loetgering, Kevin Zhou, and Roarke Horstmeyer

Pioneers

- 46 A Remembrance of David C. Joy, a True Microscopy and Microanalysis Pioneer**
 Dale E. Newbury and Cameron Varano Casasanta

Highlights from *Microscopy and Microanalysis*

- 52 Application of Synchrotron Radiation-Based Micro-Analysis on Cadmium Yellows in Pablo Picasso's *Femme***
- 52 A Three-Dimensional Reconstruction Algorithm for Scanning Transmission Electron Microscopy Data from a Single Sample Orientation**
- 53 Real-Time Integration Center of Mass (riCOM) Reconstruction for 4D STEM**

Departments

- 7 Editorial**
- 8 Carmichael's Concise Review**
- 48 Industry News**
- 50 Product News**
- 54 Picture Puzzle**

- 56 NetNotes**
- 66 Calendar**
- 69 Dear Abbe**
- 70 Index of Advertisers**



The most comprehensive source for all fields of microscopy, histology and general laboratory research... **gets an upgrade.**

It's the EMS website you've been waiting for. An all-new design with the familiar huge selection of products, services, and technical info, plus improved quoting system, shipping calculator, and customer accounts that will allow you to shop faster, stay updated on your order's status, and keep track of your previous orders.

www.emsdiasum.com

The New EMS Website...

ONLINE NOW ▶▶▶

**Electron
Microscopy
Sciences**

P.O. Box 550 • 1560 Industry Rd • Hatfield, PA 19440
Tel: (215) 412-8400 • Fax: (215) 412-8450
email: info@emsdiasum.com • www.emsdiasum.com

Featured Products

Discover High-Quality Products You Need in Your Lab

Bulk Chemicals

EMS has introduced a new line of bulk buffers and fixatives specifically for labs requiring higher volumes.

C-flat™ Holy Carbon Grids

C-flat™ is the premier holey carbon grid for cryo-transmission electron microscopy.

NIGHTSEA® Fluorescence

The NIGHTSEA® Stereo Microscope Fluorescence Adapter allows you to do any stereo microscopy for fluorescence method modification.

Access over 25,000 Products

EMS Virtual Catalogs and PDF Downloads

Full Line Catalog XIX

Our complete selection of products - your comprehensive EMS you know, now available.

Buffers and Fixatives

The leading products from EMS you know, now available.

New Product FOCUS

This catalog is frequently updated with the latest, the latest.

Cryo Supplies & Accessories

Featuring Cryogenic Film and Plates, Freeze Blocks, Cryo-

Access our PDF and Virtual Catalogs

Access Technical Information & News

Connect and Learn From the Community

View the Resource | See Our Team | Get Technical Information | Other Training Videos

Featured Content

2022 Trade Shows

October 2022: The 2022 Electron Microscopy Society of America (EMSA) meeting will be held in San Francisco, CA. The meeting will feature a variety of technical sessions, workshops, and networking opportunities.

View the list

Just In

Browse The Newest Products

EMS CryoDefender

GODEX DT200 Direct Thermal Printer

Access New Products

Our Companies

Visit All of Our Affiliate Websites

Learn More

citilife

DIATOME

EMS Microscopy Academy

NegaFile

rubis

Contract Packaging

DUMONT

idealtek

NIGHTSEA

SUMMERS OPTICAL

Access our Affiliates

