


Volume 18, Number 5

October 2012

Microscopy AND Microanalysis




CAMBRIDGE
UNIVERSITY PRESS

ISSN 1431-9276

Living up to Life

Leica
MICROSYSTEMS



NEW! Leica EM ACE Coaters



Courtesy of Mag. D Gruber, Core Facility of Cell Imaging and Ultrastructure Research, University of Vienna, Austria

Leica Microsystems Announces the Release of the **NEW** Leica EM ACE Coaters

The ultimate quality of a microscopic image comes from high-quality sample preparation. Leica Microsystems has a full range of innovative instrumentation to deliver perfect preparation for all sample materials.

NEW! Leica EM ACE. Coating systems for the preparation of samples for the Scanning Electron Microscope (SEM), Transmission Electron Microscope (TEM) and the Field Emission Scanning Electron Microscope (FE-SEM).

- **Fully Automated, One-touch Coating Systems**
- **Low Vacuum or High Vacuum with Optional Cryo**
- **Small Footprint**

Visit www.leica-microsystems.com/products/electron-microscope-sample-preparation
to see the latest sample preparation solutions!



Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

THE OFFICIAL JOURNAL OF

MICROSCOPY SOCIETY OF AMERICA
MICROANALYSIS SOCIETY
MICROSCOPICAL SOCIETY OF CANADA /
SOCIÉTÉ DE MICROSCOPIE DU CANADA
MEXICAN MICROSCOPY SOCIETY
BRAZILIAN SOCIETY FOR MICROSCOPY AND MICROANALYSIS
VENEZUELAN SOCIETY OF ELECTRON MICROSCOPY
EUROPEAN MICROBEAM ANALYSIS SOCIETY
AUSTRALIAN MICROSCOPY AND MICROANALYSIS SOCIETY
PORTUGUESE SOCIETY FOR MICROSCOPY

PUBLISHED IN AFFILIATION WITH

ROYAL MICROSCOPICAL SOCIETY
GERMAN SOCIETY FOR ELECTRON MICROSCOPY
BELGIAN SOCIETY FOR MICROSCOPY
MICROSCOPY SOCIETY OF SOUTHERN AFRICA

Editor in Chief
Editor, Biological Applications
Robert L. Price
Cell and Developmental Biology and
Anatomy
University of South Carolina
Columbia, SC 29209
e-mail: Bob.Price@uscmed.sc.edu

Editor, Materials Applications
David J. Smith
Department of Physics
School of Materials
Arizona State University
Tempe, Arizona 85287-1504
e-mail: david.smith@asu.edu

**Editor, Scanning Probe
Microscopies**
David Bell
School of Engineering & Applied
Sciences
Center for Nanoscale Systems
Harvard University
Cambridge, MA 02138
e-mail: dcb@seas.harvard.edu

Editor, Atom Probe
Thomas Kelly
Cameca Instruments, Inc.
A Business Unit of AMETEK, Inc.
Madison, WI 53711-4951
e-mail: Thomas.Kelly@ametec.com

**Editor, Light and Fluorescence
Microscopies**
Brian Herman
Cellular and Structural Biology
University of Texas at San Antonio
San Antonio, Texas 78284-7762
e-mail: hermanb@uthscsa.edu

Editor, Biological Applications
Heide Schatten
Veterinary Pathobiology
University of Missouri-Columbia
Columbia, Missouri 65211-5030
e-mail: schattenh@missouri.edu

Editor, Microanalysis
John Mansfield
Electron Microbeam Analysis Lab
North Campus, 417 SRB
University of Michigan
Ann Arbor, MI 48109-2143
e-mail: jfmjfm@umich.edu

**Editor, Correlative and Emerging
Microscopy Applications**
Vinayak P. Dravid
Materials Science and Engineering
Northwestern University
Evanston, Illinois 60208-3105
e-mail: v-dravid@northwestern.edu

Editor, Plant Biology Applications
Rosemary White
CSIRO Plant Industry
Canberra, ACT 2601, Australia
e-mail: rosemary.white@csiro.au

Special Issues and Reviews Editor
Jay Jerome
Vanderbilt University Medical Center
Nashville, TN 37232
e-mail: jay.jerome@vanderbilt.edu

Book Review Editor
Cynthia S. Goldsmith
Centers for Disease Control
Atlanta, GA 30333
e-mail: csg1@cdc.gov

M&M Program Guide Editor
Richard E. Edelmans
Miami University
Oxford, OH 45056
e-mail: edelmare@muohio.edu

Proceedings Editor
John Shields
University of Georgia
Athens, GA 30602
e-mail: jpshield@uga.edu

Editorial Board

Ralph Albrecht	<i>University of Wisconsin, Madison, Wisconsin</i>
Barry Carter	<i>University of Connecticut, Storrs, Connecticut</i>
Wah Chiu	<i>Baylor College of Medicine, Houston, Texas</i>
Niels de Jonge	<i>INM Institute for New Materials, Saarbrücken, Germany</i>
Alberto Diaspro	<i>University of Genoa, Italy</i>
Elizabeth Dickey	<i>Pennsylvania State University, University Park, Pennsylvania</i>
Alwyn Eades	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Mark Ellisman	<i>University of California at San Diego, San Diego, California</i>
Pratibha Gai	<i>University of York, United Kingdom</i>
Marija Gajdardziska-Josifovska	<i>University of Wisconsin-Milwaukee, Milwaukee, Wisconsin</i>
Dale Johnson	<i>University of South Florida, Tampa, Florida</i>
Paul Kotula	<i>Sandia National Labs, Albuquerque, New Mexico</i>
William Landis	<i>University of Akron, Akron, Ohio</i>
Eric Lifshin	<i>SUNY at Albany, Albany, New York</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Dale Newbury	<i>National Institute of Standards and Technology, Gaithersburg, Maryland</i>
Jean-Paul Revel	<i>California Institute of Technology, Pasadena, California</i>
Conly Rieder	<i>Wadsworth Center, Albany, New York</i>
John Silcox	<i>Cornell University, Ithaca, New York</i>
Nestor Zaluzec	<i>Argonne National Laboratory, Argonne, Illinois</i>

Editorial Board Representatives from Affiliated Societies

Ian Anderson	<i>NIST, Gaithersburg, Maryland (MAS)</i>
Gautam Kumar Dey	<i>Bhabha Atomic Research Centre (EMSI)</i>
Gema Gonzalez	<i>Venezuelan Institute for Scientific Investigation (Venezuela)</i>
Michael Robertson	<i>Acadia University, Wolfville, Nova Scotia (Canada)</i>
Brendan Griffin	<i>University of Western Australia (AMMS)</i>
Guillermo Solorzano	<i>Pontificia Universidade Catolica, Rio de Janeiro (Brazil)</i>
Clive Walker	<i>Institute for Transuranium Elements, Karlsruhe (EMAS)</i>
Miguel Yacamán	<i>Mexico Institute for Nuclear Research (Mexico)</i>
Henrique Almeida	<i>Universidade do Porto (Portugal)</i>

Founding Editor

Jean-Paul Revel	<i>California Institute of Technology, Pasadena, California</i>
-----------------	---

Previous Editors-in-Chief

Dale Johnson	<i>University of South Florida, Tampa, Florida</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>

This journal is part of the **Cambridge Journals Online** service. Access to online tables of contents and article abstracts is available to all researchers at no cost. Access to full-text articles online is provided to those with online subscription. Online subscriptions must be activated. Once your subscription is activated, free access to past, present, and forthcoming articles is available at:

Microscopy and Microanalysis website: journals.cambridge.org/MAM.

Instructions for authors submitting manuscripts may be found at journals.cambridge.org/MAM. Select "Further Information" then select "Instructions for Contributors." An abbreviated version of these instructions will be published in the first issue (February) of each volume.



Cooling Stages
Recirculating Heaters and Chillers
Sputter Coaters
SEM/TEM Carbon Coaters
Vacuum Evaporators
Glow Discharge Systems
RF Plasma Etchers/
Plasma Reactors
Critical Point Dryers
Freeze Dryers
Cryo-SEM Preparation Systems
Vacuum Pumps & Accessories
Evaporation Supplies
and more...

well equipped...

Electron Microscopy Sciences is pleased to announce our new full line catalog, your comprehensive source for high-end Vacuum Equipment. EMS is committed to providing the highest quality products along with competitive pricing, prompt delivery and outstanding customer service.



and more...

Not just Vacuum Equipment, EMS also offers:

Laboratory Microwave Ovens • Automated Tissue Processors • Oscillating Tissue Slicers • Vibrating Microtomes • Rapid Immersion Freezers • Tissue Choppers • Desiccators and Desiccants • Centrifuges, Tubes, and Racks • Stirrers, Stirring Hotplates, and Digital Hotplates • Stirring Bars, Stirring Rods, and Hand Mixers • Vortex Mixers, Microplate Mixers, and Magnetic Stirrers • Tissue Rotators, Mixer Vortexes, and Rotator/Rockers • Dri Baths • Oven/Incubators • Cooling Chambers • Ultraviolet Lamps • Lab Jacks


For catalog requests, please visit our website at www.emsdiasum.com



**Electron
Microscopy
Sciences**

Electron Microscopy Sciences
P.O. Box 550 • 1560 Industry Rd. • Hatfield, Pa 19440
Tel: (215) 412-8400 • Fax: (215) 412-8450
email: sgkcock@aol.com • www.emsdiasum.com

The SPI Supplies Family of Instruments



Excellence in sample preparation
just a click away... 2spi.com

Your results will never be better than your sample preparation. See how SPI Supplies can help you deliver the highest quality results for all your SEM/EDS, TEM and FESEM applications.

1. Osmium Plasma Coaters for FESEM Applications
2. SPI-MODULE™ Sputter/Carbon Coater Module
3. Gentle Mill™ Ion Milling System
4. Plasma Prep™ III Solid State Plasma Cleaner for cleaning TEM holders

5. Plasma Prep™ X Parallel Plate Plasma Etcher
6. Plasma Prep™ III Plasma Etcher with PPIII Process Controller
7. Vacu Prep™ II Turbo Pump Evaporation System
8. SPI-DRY™ Critical Point Dryer
9. Precision Spin Coater Spin coater



SPI Supplies Division of **STRUCTURE PROBE, Inc.**






P.O. Box 656 • West Chester, PA 19381-0656 USA
Phone: 1-610-436-5400 • 1-800-2424-SPI (USA and Canada) • Fax: 1-610-436-5755 • 2spi.com • E-mail: sales@2spi.com





High Definition Digital
TEM Cameras with
1 to 16 Megapixels

- AMT SOLUTIONS
- Life Science Cameras
- Material Science Cameras
- Easy To Use Software
- Reliability and Services
- TEM Integration
- Extensive Support




3 Electronics Avenue, Danvers, MA 01923 • Phone: 978.774.5550 • www.amtimaging.com



Explore the Unseen

FEI's annual owner's image contest
In partnership with National Geographic



FEI is pleased to partner with National Geographic on a giant-screen film tentatively titled "Invisible Worlds". At the same time, we're excited to announce this year's FEI Image Contest, "Explore the Unseen". Inspired by the upcoming film, the FEI Image Contest offers owners and users an opportunity to explore their creativity and share their images with National Geographic's worldwide audience.

For more information go to www.fei.com/2012-image-contest/

All images courtesy of FEI and can be found at www.flickr.com/photos/fei_company

Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

Volume 18, Number 5
October 2012

Review Article

What Remains to Be Done to Allow Quantitative X-Ray Microanalysis Performed with EDS to Become a True Characterization Technique? 915

Raynald Gauvin

TECHNIQUES AND EQUIPMENT DEVELOPMENT

Refining Spatial Distribution Maps for Atom Probe Tomography via Data Dimensionality Reduction Methods 941

Santosh K. Suram and Krishna Rajan

On the Use of Simulated Field-Evaporated Specimen Apex Shapes in Atom Probe Tomography Data Reconstruction 953

David J. Larson, Brian P. Geiser, Ty J. Prosa, and Thomas F. Kelly

Atom Probe Tomography Simulations and Density Functional Theory Calculations of Bonding Energies in Cu₃Au 964

Torben Boll, Zhi-Yong Zhu, Talaat Al-Kassab, and Udo Schwingenschlögl

Atom Probe Tomography of Nickel-Based Superalloys with Green or Ultraviolet Lasers: A Comparative Study 971

Yaron Amouyal and David N. Seidman

Atomic Resolution Phase Contrast Imaging and In-Line Holography Using Variable Voltage and Dose Rate 982

Bastian Barton, Bin Jiang, ChengYu Song, Petra Specht, Hector Calderon, and Christian Kisielowski

Strain Measurement in Semiconductor Heterostructures by Scanning Transmission Electron Microscopy 995

Knut Müller, Andreas Rosenauer, Marco Schowalter, Josef Zweck, Rafael Fritz, and Kerstin Volz

Improving Fabrication and Application of Zach Phase Plates for Phase-Contrast Transmission Electron Microscopy 1010

Simon Hettler, Björn Gamm, Manuel Dries, Nicole Frindt, Rasmus R. Schröder, and Dagmar Gerthsen

Fast Amplitude Modulation Detector for Scanning Force Microscopy with High Q Factor 1016


Yizhi Shi and Qingyou Lu

Partial Acceptor Photobleaching-Based Quantitative FRET Method Completely Overcoming Emission Spectral Crosstalks 1021

Huali Li, Huaina Yu, and Tongsheng Chen

Visualization of ATP with Luciferin-Luciferase Reaction in Mouse Skeletal Muscles Using an "In Vivo Cryotechnique" 1030

N. Terada, Y. Saitoh, S. Saitoh, N. Ohno, K. Fujishita, S. Koizumi, and S. Ohno



On the Cover: Synchrotron tomography of fossil bone. For further information see Sanchez et al. pp. 1095–1105.

Microscopy and Microanalysis website: <http://www.journals.cambridge.org/MAM>
Indexed in Chemical Abstracts, Current Contents, BIOSIS, and MEDLINE (PubMed)

Electron Tomography of HEK293T Cells Using Scanning Electron
Microscope–Based Scanning Transmission Electron Microscopy 1037
*Yun-Wen You, Hsun-Yun Chang, Hua-Yang Liao, Wei-Lun Kao, Guo-Ji Yen,
Chi-Jen Chang, Meng-Hung Tsai, and Jing-Jong Shyue*

Sample Preparation Induced Artifacts in Cryo-Electron Tomographs 1043
*Pavel Plevka, Anthony J. Battisti, Dennis C. Winkler, Kaspars Tars,
Heather A. Holdaway, Carol M. Bator, and Michael G. Rossmann*

DNA Base Identification by Electron Microscopy 1049
*David C. Bell, W. Kelley Thomas, Katelyn M. Murtagh, Cheryl A. Dionne,
Adam C. Graham, Jobriah E. Anderson, and William R. Glover*

BIOLOGICAL APPLICATIONS

Aging Increases Nuclear Chromatin Entropy of Erythroid Precursor Cells
in Mice Spleen Hematopoietic Tissue 1054
Igor Pantic, Senka Pantic, and Jovana Paunovic

Consequences of a Fat Diet in the Distribution of Minerals within Pancreatic
Tissues of Rats and Rabbits 1060
M. Dolores Ynsa, Ren Minquin, Reshmi Rajendran, Teresa Pinheiro, and Frank Watt

Ultrastructural and Biochemical Alterations in Rats Exposed to Crude Extract
of *Carex baccans* and *Potentilla fulgens* 1067
Bishnupada Roy, Bikash Ranjan Giri, Mitali Chetia, and Ananta Swargiary

Microscopy and Microanalysis of Hematological Parameters in Common Carp,
Cyprinus carpio, Inhabiting a Polluted Lake in North East India 1077
Bashida Massar, Sudip Dey, R. Barua, and K. Dutta

Detection of Biofilm-Grown *Aspergillus fumigatus* by Means of Atomic Force
Spectroscopy: Ultrastructural Effects of Alginate Lyase 1088
*M. Papi, A. Maiorana, F. Bugli, R. Torelli, B. Posteraro, G. Maulucci, M. De Spirito,
and M. Sanguinetti*

Three-Dimensional Synchrotron Virtual Paleohistology: A New Insight into
the World of Fossil Bone Microstructures 1095
*Sophie Sanchez, Per E. Ahlberg, Katherine M. Trinajstic, Alessandro Mirone,
and Paul Tafforeau*

Three-Dimensional Morphology of Microdamage in Peri-Screw Bone:
A Scanning Electron Microscopy of Methylmethacrylate Cast Replica 1106
Lei Wang, Jin Shao, Tingjun Ye, Lianfu Deng, and Shijing Qiu

Micro Energy-Dispersive X-Ray Fluorescence Mapping of Enamel and
Dental Materials after Chemical Erosion 1112
*Luís Eduardo Silva Soares, Rodrigo de Oliveira, Sídnei Nahórny,
Ana Maria do Espírito Santo, and Airton Abrahão Martin*

MATERIALS APPLICATIONS

Three-Dimensional Chemistry of Multiphase Nanomaterials by Energy-Filtered
Transmission Electron Microscopy Tomography 1118
*Lucian Roiban, Loïc Sorbier, Christophe Pichon, Pascale Bayle-Guillemaud,
Jacques Werckmann, Marc Drillon, and Ovidiu Ersen*

Nanometer Scale Tomographic Investigation of Fine Scale Precipitates in a
CuFeNi Granular System by Three-Dimensional Field Ion Microscopy 1129
*Sophie Cazottes, François Vurpillot, Abdelsem Fnidiki, Dany Lemarchand,
Marcello Baricco, and Frederic Danoix*


Disassembling Glancing Angle Deposited Films for High-Throughput, Single-Post
Growth Scaling Measurements 1135
*Joshua Morgan Arthur Siewert, Joshua Michael LaForge, Michael Thomas Taschuk,
and Michael Julian Brett*

Contents continued

Optoelectronic Properties of InAlN/GaN Distributed Bragg Reflector Heterostructure Examined by Valence Electron Energy Loss Spectroscopy <i>A. Eljarrat, S. Estradé, Ž. Gačević, S. Fernández-Garrido, E. Callejo, C. Magén, and F. Peiró</i>	1143
Direct Observations of Confined Layer Slip in Cu/Nb Multilayers <i>Nan Li, Jian Wang, Amit Misra, and Jian Yu Huang</i>	1155
The Use of μ CT and ESEM in the Study of the Osmosis-Induced Water Uptake by Eurobitum Bituminized Radioactive Waste <i>An Mariën, Elie Valcke, Nele Bleyen, Maarten Van Geet, and Martine Wevers</i>	1163
Microstructural Characterization of Consolidant Products for Historical Renders: An Innovative Nanostructured Lime Dispersion and a More Traditional Ethyl Silicate Limewater Solution <i>Giovanni Borsoi, Martha Tavares, Rosário Veiga, and Antonio Santos Silva</i>	1181
Microstructural Analysis of Iberian Expanded Clay Aggregates <i>J. Alexandre Bogas, António Mauricio, and M.F.C. Pereira</i>	1190

Count on Us!

Best EDS Performance with the **NEW** Slim-line XFlash® 6




You can count on the **NEW** XFlash® SDD generation:

- Best solid angle – optimum geometry and active areas from 10 mm² to 100 mm²
- Best throughput – up to 600,000 cps output at 1,500,000 cps input
- Best energy resolution – 121 eV at Mn K α , 47 eV at F K α , 38 eV at C K α (FWHM, exceeds ISO 15632:2002 requirements)

www.bruker.com

Preparation Equipment and Microscopy Supplies

The single source for All your microscopy supplies and specimen preparation equipment.



- Vacuum Coating Systems
- Calibration Standards
- PELCO® easiGlow™ Glow Discharge Unit
- SEM Sample Holders and Mounts
- Silicon Nitride TEM Membranes
- PELCO BioWave Pro® Tissue Processor
- TEM Support Films



- AFM Supplies
- Quality Laboratory Tweezers
- Vacuum Pick-up Systems
- Digital Stereo Microscopes
- Conductive Adhesives
- FIB Supplies



Complete line of compact Cressington EM Sample Coaters.

 **TED PELLA, INC.**
Microscopy Products for Science and Industry

sales@tedpella.com

800-237-3526

www.tedpella.com


Unrivaled EDS for TEM

0.98sr 100mm²

- Exponentially enhances elemental mapping for nano-area analysis TEM
- Automatic retractable design




CENTURIO
Large Solid Angle SDD for TEM



*High sensitivity
for fast mapping
at atomic resolution*

*Seamless chemical mapping and
data collection for S/TEM-TEM-
SEM-EDS. (256 x 256 pixels.
Total acquisition time: 1 min. 13 sec.)*



*View our real-time
DRAM analysis video
at www.jeolusa.com/UnrivaledEDS*

JEM-2800



JEM-ARM200F



JEOL

Another
Extreme Imaging
Solution

Global Solutions Provider for Advanced Technology
www.jeolusa.com • salesinfo@jeol.com
978-535-5900

Find us on Facebook and Twitter @jeolusa

NANOANALYSIS S/TEM